**SICUREZZA**  
 Conforme a:  
 EN60601-1  
 EN60601-1-2  
 EN60601-2-2

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**TWIN OPERATIVE ACTIVITY**

- Two output channels for monopolar cutting
- Two independent setting and output channels for the most frequent monopolar coagulation

**CONTROLLED CUTTING PROCEDURES**

- Enhanced cutting current for eschar free incision
- Cutting current with selectable coagulating property for bleeding controlled dissection
- Low capacitive coupling by Slow Blend current for endoscopic procedures

**TWIN CONTROLLED COAGULATIVE PROCEDURES**

- Speedy\* coagulation and Deep\*\* coagulation available to two operators
- Spray coagulation available by hand switch or foot switch

**WIDE BIPOLAR FEATURES**

- Bipolar cutting current with selectable bleeding control
- Slow Blend cutting current
- Separately activable Automatic Start and Automatic Stop coagulation

**SURGEON AND PATIENT TOTAL SAFETY**

- Acceptability of contact impedance between patient body and return electrode continuously watched by Skin Plate Electronic Control (SPEC)
- Excessive power delivery time warned
- H.F. Leakage currents continuously monitored
- Mains voltage variations automatically compensated before to start the delivery of HF power

**PERSONALIZED PROCEDURES**

- Independent choice of coagulative currents and power levels by each of two operators
- Ten personal configurable and immediately recallable procedures

**EXTENTION DRIVER AND INTERLOCK SIGNALS**

- Drived output connection to external Argon Plasma unit and Smoke Evacuator unit

**AUTOMATIC SELF DIAGNOSIS**

- Automatic diagnostic control of the actual complete functionality of the electronic unit and the connected accessories and coded report of found failures

\*Speedy coagulation also said Fulgurate or Forced coagulation  
 \*\*Deep coagulation als said Pin Point, Dessiccate or Soft coagulation

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**STANDARD COMPOSITION**

Power supply cable 5 m. SIEMENS-IEC (1 pcs)
Cable for connection neutral plate (1 pcs)
Disposable split neutral electrode (5pcs)
Blade electrode 7 cm (3 pcs)
Blade electrode 16 cm (3 pcs)
Needle electrode 7 cm (3 pcs)
Ball electrode 6 cm (3 pcs)
Reusable handpiece with finger switches (1 pcs)
Reusable handle with switches (1 pcs)
Double water-proof foot-switch (2 pcs)
Manual of instructions
Electrode cleaning sponge (1 pcs)

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HIGH POWER  
 MONOPOLAR AND BIPOLAR  
 ELECTROSURGICAL UNIT

CHIRURGIA AMBULATORIALE  
 CHIRURGIA GENERALE  
 CHIRURGIA PEDIATRICA  
 CHIRURGIA PLASTICA  
 CHIRURGIA TORACICA  
 CHIRURGIA VASCOLARE  
 DERMATOLOGIA  
 ENDOSCOPIA  
 GASTROENTEROLOGIA  
 GINECOLOGIA  
 NEUROCHIRURGIA  
 ORTOPEDIA  
 OTORINOLARINGOIATRIA  
 PNEUMOLOGIA  
 UROLOGIA  
 VETERINARIA

CASUALTY  
 GENERAL SURGERY  
 PEDIATRIC SURGERY  
 PLASTIC SURGERY  
 THORAX SURGERY  
 VASCULAR SURGERY  
 DERMATOLOGY  
 ENDOSCOPY  
 GASTROENTEROLOGY  
 GYNECOLOGY  
 NEUROSURGERY  
 ORTHOPEDICS  
 OTORHINOLARYNGOLOGY  
 PNEUMOLOGY  
 UROLOGY  
 VETERINARY

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Adjustment of the levels of power of the functions of the monopolar section CUT:  
**CUT**, for pure cut or with coagulative properties, suited for dissection;  
**ENHANCED**, for incision with reduced production of eschar;  
**BLEND**, for endoscopic procedures minimally invasive, with low capacity coupling through the surrounding tissues  
Selection of the function to deliver.  
Visualization of the level of power of the selected function, that can be activated through command of the connected accessories.

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
Observation of the state of mains power supply and electronic interrupters for lighting and turning off of the unit. The fluctuations of the main power supply voltage are automatically compensated before to start the delivery of the HF current.

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User interface for:  
-Planning of the time of delay of the beginning of the automatic bipolar coagulation. This time can be programmed to 0,5s, 1s or 1,5s.  
-Reading of last 101 errors eventually happened.  
-Recording of 10 personalized procedures; selection and recall of the same. In the personalized procedures it is possible to record both the selected monopolar and bipolar functions and the relative levels of power.

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Connection of normal or split neutral electrode. Indication and acceptance of the plate to patient impedance. The circuit gives a signal of alarm and inhibits the delivery of currents if the safety impedance limits have exceeded.

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Adjustment of the levels of power of the functions of the monopolar section COAG1:  
**SPEEDY**, for rapid coagulation;  
**DEEP**, for coagulation and deep desiccation.  
Selection of the function to deliver.  
Visualization of the level of power of the selected function, that can be activated through command of the connected accessory.

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
Adjustment of the levels of power of the functions of the monopolar section COAG2:  
**SPEEDY**, for rapid coagulation;  
**DEEP**, for coagulation and deep desiccation;  
**SPRAY**, for wide superficial coagulation without contact;  
**ARGON SPRAY**, supported by unity Argon-Plasma.  
Selection of the function to deliver.  
Visualization of the level of power of the selected function, that can be activated through command of the connected accessories.

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Adjustment of the levels of power of the bipolar CUT:  
**CUT**, for pure cut or with coagulative property, suited for dissection;  
**BLEND**, low voltage current, proper to endoscopic operations, thanks to the low risk of power transfer by capacitive coupling to the surrounding tissues.  
Digital visualization of the bipolar cut output power.

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Adjustment of the power level of bipolar coagulation:  
Choice of the modality of delivery of the current can be also selected with automatic beginning (START) and/or automatic end (STOP). The automatic beginning of the delivery of current can take place 0,5s, 1s or 1,5s after the bipolar forceps contact with the patient, in order to improve and to favor the process of coagulation.  
Visualization of the power that can be activated by foot-switch.

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Output for the connection of the bipolar accessories with command of activation through double pedal for the selection of the functions.

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The bright signal OVT points out exceeding of the maximum time, recommended for delivery of current, to avoid excessive heating of the treated area.

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Two output channels, reciprocally locked, with prevalence of the first activated. The connection of two handpieces with switch button and of a double pedal allows the personalized selection of the functions.  
Channel CUT/COAG1 can deliver, through commands from handpiece, the cut currents or the coagulation currents available in COAG1.  
Channel CUT/COAG2 can deliver, with commands by handpiece or by foot-switch, the cut currents or the coagulation currents available in COAG2.  
The bright signal OVT points out exceeding of the maximum time, recommended for delivery of current, to avoid excessive heating of the treated area.

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DP258/GBb MADE IN ITALY © 2006 LED SpA
**LED SpA**