





MAHE-Medizintechnik Made in Germany



Urethrotom
Urethrotom



Urethrotom (Urologie) Urethrotome (Urology)




Optiken Telescopes

Blickrichtung
Direction of view

0°	35.001.00
30°	35.001.30
70°	35.001.70

AUTOCLAVE
134° C / 273° F




Arbeitseinheit Resektoskop (passiv) Working unit (passive)

36.005.01


HF-Anschlußkabel für Arbeitseinheiten HF-connecting cord for working units

36.005.10




Lichtleit-Kabel Glasfaser Cold light cable glas fiber

Ø	Nutzlänge Working length	
3,5 mm	180 cm	11.000.01
3,5 mm	230 cm	11.000.02
4,8 mm	180 cm	11.000.03



Urethrotom-Schaft Urethrotome sheath


37.006.01



Obturator allein Obturator only

36.006.03

Projektor-Adapter




Storz




11.001.00

Olympus




11.001.02

Wolf/Dyonics




11.001.04

ACM/British standard




11.001.06



Führungsrohr, seitlich offen, zum Legen eines Ballonkatheters Guide tube, half-round, to insert a ballon catheter

37.007.01

Endoskop-Adapter



Storz



11.001.20

Wolf



11.001.22

ACM/British standard




11.001.24

Stryker




11.001.26



Schaft für kontinuierliche Spülung / Absaugung Sheath for continuous suction / irrigation

37.007.05




Strikturmesser gerade
Cold knives straight

37.015.01

Lagerungs- und Sterilisierbehälter
Protection and sterilization tray




37.015.00



Strikturmesser halbrund
Cold knives half round

37.015.02

Andere Ausführungen auf Anfrage
 Other models on request



Ballon 15 ml, Luer-Ventil
15 ml balloon, Luer valve

CH/Fr	
12	AA3812
14	AA2814
16	AA3816
18	AA3818
20	AA3820

Uretherkatheter, zylindrisch
Cylindrical ureteral catheter

CH/Fr	
3	AC5303
4	AC5304
5	AC5305
6	AC5306
7	AC5307
8	AC5308
9	AC5309
10	AC5310
11	AC5311

CH/Fr	
12	AA3712
14	AA2714
16	AA3716
18	AA3718
20	AA3720
22	AA3722



MAHE-Medizintechnik
Föhrenstrasse 45 D-78532 Tuttlingen/Germany
Tel. (0 74 61) 1 31 47 Fax (0 74 61) 7 98 23